

## CONVENTIONAL, ISO COLOUR-CODED, ISO TAPER, AND MULTI-TAPER PAPER POINTS.



### Application:

- For drying root canals so that the sealing and filling products will adhere better to the canal walls.

## Why choose R&S paper points?

- + Highly absorbent**  
Perfect quick drying of canals.  
Won't tear even when wet.
- + Uniformity**  
- Taper and uniform consistency.  
- Smooth, rounded points.  
- Perfect for preparing canals.
- + Production quality**  
All of R&S's points are in compliance with and surpass Standard ISO 7551 and ADA specifications, and they are individually inspected by laser.
- + Wide selection**  
A wide variety of sizes to meet all of your needs.  
Packaged in sterile leakproof packaging.



## Advice from our dental experts

- After drying the canal, keep a paper point in place before filling.

## Technical specifications and instructions for use

ISO STANDARD	7551 - DENTAL ABSORBENT POINTS			
Type	CONVENTIONAL	ISO COLOUR-CODED	ISO TAPER	ISO MULTI-TAPER
Use	With all types of sequences	With all conventional instruments	With all instruments with constant tapers, such as PROFILE®, CMA®, HERO 642®, HERO SHAPER® and R-endo®	With all instruments with variable tapers, such as PRO Taper®, Go Taper Universal & Flex®
Shelf life (in years)			4	
Taper (%)	Variable	2	4 and 6	Variable
ISO (Ø)	 XX-Fine X-Fine Fine Medium Large X-Large Assorted	 15 20 25 30 35 40 45 50 55 60 Assorted 15-40 Assorted 45-80	 15 20 25 30 35 40 45 Assorted 15-40 Assorted 45-80	 F1 F2 F3 Assorted F1-F3 Assorted F4-F5
Length (mm)		28		30
Sterilisation		By gamma rays		
Packaging		Box of 200 points		Box of 100 points

## Additional products


Dams with or without latex




Screw-posts



Nerve broaches



Easyaccess Niti



K-files, K-reamers and H-files



Endo EZ burs



Paper points : class IIa medical devices. CE0197. Certifying body: TUV RHEINLAND. Manufactured by DIADENT GROUP INTERNATIONAL. Dams : class I medical devices. CE marking. Manufactured by FOUR D RUBBERCO LTD. Nerve broaches, K-files, K-reamers, H-files : class I medical devices. CE marking. Manufactured by CFPM. Screw posts : class IIa medical devices. CE1639. Certifying body: SGS BELGIUM. Manufactured by NORDIN HARALD. Endo EZ burs: Class IIA medical devices. CE0123. Certifying body : TUV PRODUCT SERVICE. Manufactured by PRIMA DENTAL GROUP.

Images are for illustrative purposes only.  
 GACD - 25 rue Bleue, 75009 PARIS - France 435 051 776 R.C.S PARIS - SASU.  
 Tel.: +33 (0)1 48 01 27 14 - Fax: +33 (0)1 45 23 51 89.  
 export@rs-dentalproducts.com - www.rs-dentalproducts.com.

